

Code; RP802-xxx

Lot; xxxxx

Size; xxxxxunits



ProteinExpress

PNGaseF (N-Glycosidase F)

Supplied Reagents

- PNGaseF
- 10 X PNGaseF Reaction Buffer
- 10 X Denaturing Solution
- 10% NP-40

Concentration : 500,000units/mL

Storage : -20 °C

Description : PNGaseF(N-GlycosidaseF) is a recombinant glycosidase from *Elizabethkingia miricola* (*Flavobacterium Meningosepticum*).

Storage Buffer :

20 mM Tris-HCl (pH7.5 at 25°C)
50 mM NaCl
5 mM EDTA
50% Glycerol

10 X PNGaseF Reaction Buffer :

0.5 M Sodium Phosphate(pH7.5)

10 X Denaturing Buffer :

5% SDS, 0.4M DTT

Source : Recombinant protein with a N-terminal DYKDDDDK tag and 6xHis tag, expressed in *E.coli*.

Molecular Weight : 36kDa.

Purity : > 95%, as determined by SDS-PAGE visualized by CBB stain.

Unit Definition : One unit is defined as the amount of enzyme required to remove >95% of the carbohydrate from 10µg of denatured RNase B in 1hour at 37°C in a total reaction volume of 10µl.

Standard Procedure :

1. Combine 1-20 µg of glycoprotein, 1 µl of 10X Denaturing Buffer and H₂O (if necessary) to make a 10 µl total reaction volume.
2. Denature glycoprotein by heating reaction at 100° C for 10 minutes.
3. Make a total reaction volume of 20 µl by adding 2 µl of 10X PNGase Reaction Buffer, 2 µl of 10% NP-40, H₂O and 1-5 µl PNGaseF.
4. Incubate at 37 °C for 1 hour

References :

- 1)Maley, F. et al. (1989). *Anal. Biochem.* 180, 195-204.
- 2)Tretter, V. et al. (1991). *Eur. J. Biochem.* 199, 647-652.
- 3)Plummer, T.H. Jr. and Tarentino, A.L. (1991). *Glycobiology.* 1, 257-263.

For Research Use Only. Not for use in diagnostic procedures

ProteinExpress Co., Ltd.

2-1-6, Kazusakamatari, Kisarazu-shi, Chiba 292-0818, Japan

Tel: +81-438-52-0112, Fax: + 81-438-52-0113

E-mail; service@proteinexpress.co.jp

URL; <http://www.proteinexpress.co.jp>